

Genotyping sgRNA flies with pCFD4 and pCFD3 inserts

Prepared by Rong Tao and Jonathan Zirin

1. Crude genomic DNA preparation:

SB buffer:

10mM Tris-HCl pH8

1mM EDTA

25mM NaCl

Proteinase K:

Stock conc. (100X): 20mg/ml

Working conc.: 200ug/ml (actually 20ug/ml works fine)

add Proteinase K to SB buffer just before experiment

- a. place a single fly in an Eppendorf tube
- b. add 50ul SB buffer
- c. homogenize
- d. incubate at 37°C · 30min
- e. spin: 16,000g · 2 min
- f. transfer the supernatant to PCR tube
- g. incubate 95°C · 3min
- h. store at 4°C (for longer term, store in -20°C freezer)

2. Genotyping PCR:

Use these primers for pCFD4 to amplify a band of 731bp

pCFD4-F: 5'-GACACAGCGCGTACGTCCTTCG-3'

pCFD4-R: 5'-ACTCTCAGGCTCCAGGTAGG-3'

Use these primers for pCFD3 to amplify a band of ~500bp

pCFD3-F: 5'-ACGTTTTATAACTTATGCCCCCTAAG-3'

pCFD3-R: 5'-GCCGAGCACAATTGTCTAGAATGC-3'

- a. PCR rxn

2X GoTaq Green	15 ul
H ₂ O	11.25
FW primer (10uM)	0.625
RV primer (10uM)	0.625
Genomic DNA	2.5 ul
- b. PCR program for pCFD3

Step1:	95°C, 2 min
Step2(35 x):	95°C, 30 sec
	50°C, 30 sec

72°C, 1 min
Step3: 72°C, 10 min
Step4: 4°C, ∞

c. PCR program for pCFD4

Step1: 95°C, 2 min
Step2(35 x): 95°C, 30 sec
55°C, 30 sec
72°C, 1 min
Step3: 72°C, 10 min
Step4: 4°C, ∞

3. Run gel:

- a. run 4 ul PCR product on agarose gel
- b. check if PCR produces the right size of band and enough DNA for doing PCR purification

4. PCR purification:

- a. use QIAquick PCR Purification Kit
- b. elute in ~50 ul H₂O
- c. purified DNA should be 10~20ng/ul

5. Sequence purified PCR product:

- a. For pCFD3 use pCFD3-F primer
- b. For pCFD4 use pCFD4-F primer

6. Sequencing data analysis:

- a. for pCFD3 the sequence is: GTCG + sgRNA + GTTTTAGAGC
- b. for pCFD4 the sequence is: AACTTC + sgRNA1 + 496bp + sgRNA2 + GTTTT